



People.Health.Care.

ReadyWrap® Fusion Kit **NEW**

A clinically-effective wrap and compressive liner — together at last.



- Kit includes:
- ☐ ReadyWrap® Calf
 - ☐ ReadyWrap® Fusion Liner


ReadyWrap® Fusion Kit

Clinically-effective and easy-to-use.

ReadyWrap® Fusion Kit consists of an adjustable compression wrap plus liner with compressive foot and ankle for transitioning out of therapy into self-management during the daytime, nighttime, and during high or low activity.


 Off-the-Shelf

 Day or Night

 High or Low Activity

Indications

- Lymphedema
- Edema
- Leg ulcers
- Venous insufficiency
- Leg fatigue and/or tiredness



Low-stretch material

Mimics short-stretch bandaging — offering high working and low resting compression for day or night use.



Color-coded straps

Provides easy step-by-step donning guidance.



50% Overlapping straps

Helps prevent gapping.¹



Blocked support spine

Helps maintain garment form to support all-day wear.¹




Continuous compression

Garments work together to deliver therapeutic support from foot to just below knee.



Targeted compression

Delivers 8-15mmHg in the foot and ankle that easily fits in a shoe.



ReadyWrap Fusion Kit **NEW**

kit includes ReadyWrap Calf garment in Beige or Black and single Black ReadyWrap Fusion Liner


Size	Circumference			Length		Wrap Color	REF	HCPCS			
	A ^{1c}	B ^c	C ^c	AY ^L	BD ^L						
Regular											
S	21-24	18-23	28-38	< 25	up to 30	Beige	154894	A6545 [†]			
						Black	154895				
M	25-28	20-28	33-43			Beige	154896				
						Black	154897				
L	29-32	25-33	38-48			Beige	154898				
						Black	154899				
XL	33-36	30-38	43-58			Beige	154900				
						Black	154901				
XXL	37-40	35-43	53-68			Beige	154902				
						Black	154903				
Long											
S	21-24	18-23	28-38			≥ 25	up to 35		Beige	154904	A6545 [†]
				Black	154905						
M	25-28	20-28	33-43	Beige	154906						
				Black	154907						
L	29-32	25-33	38-48	Beige	154908						
				Black	154909						
XL	33-36	30-38	43-58	Beige	154910						
				Black	154911						
XXL	37-40	35-43	53-68	Beige	154912						
				Black	154913						

ReadyWrap Fusion Liners **NEW**

Below Knee liner with focused 8-15 mmHg in foot and ankle, sold in pairs, available in Black only

Size	Circumference			Length		REF
	A ^{1c}	B ^c	C ^c	AY ^L	YD ^L	
Regular						
S	21-24	18-23	28-38	< 25	< 40	153523
M	25-28	20-28	33-43			153524
L	29-32	25-33	38-48			153525
XL	33-36	30-38	43-58			153526
XXL	37-40	35-43	53-68			153527
Long						
S	21-24	18-23	28-38	≥ 25	≥ 40	153528
M	25-28	20-28	33-43			153529
L	29-32	25-33	38-48			153530
XL	33-36	30-38	43-58			153531
XXL	37-40	35-43	53-68			153532

Measurement reference



Medicare Reimbursable

ReadyWrap Fusion Kit is covered under code **A6545[†]**.

- 6-month warranty
- Hand wash, lay flat to dry



People.Health.Care.

How to Order

1 Measure

Take the circumference and length measurements using the sizing chart and diagram provided.

2 Order

Available through your preferred garment supplier:

- Local or National DME
- National provider
- Online distributor

 www.lohmann-rauscher.us/where-to-buy

Dealers can purchase direct from L&R USA via phone, fax, or online store.

Material content

Calf: Nylon (polyamide), polyurethane, spandex (elastane)
Not made with natural rubber latex.



Liner: 95% nylon (polyamide), 5% spandex (elastane)
Not made with natural rubber latex.



References:

1. Ehmann, S., and Bock, K. (Jan. 2018), "Case report to demonstrate the need for selection criteria for optimal adjustable Velcro wrap prescription", The Journal of Advanced Wound Care, North American Supplement, vol. 27, no. 1

[†] For presence of open wound. Refer to LCD requirement for more information.

